



# In house verification of piston pipettes

Covid-19



# In House Verification of Piston Pipettes

## A handy Metrological Approach in the COVID-19 Crisis

### Authors

*Dr Zoe Metaxiotou (EIM / Greece)*

*Elsa Batista (IPQ / Portugal)*





# Introduction



The unprecedented conditions public health institutions experience due to COVID-19 pandemic crisis have forced hospital laboratories to ***perform on a daily basis a vast amount of diagnostic tests using pipettes for the handling of the samples.***





## COVID-19 situation


**Correct measurement of sample volume depends on the calibration status of the pipette**



**In pandemic situations the analytical laboratories may not have the time to send the pipettes to a calibration laboratory**



# COVID-19 situation








# COVID-19 situation

This document intends to provide hospital laboratories with a brief, handy “**Rescue Manual**” for performing a quick **in-house verification** of the measuring status of pipettes in order to identify the ones which do not perform safely within specifications





# Calibration<sup>1</sup> of piston operated pipettes

## Routine calibration at 3 points:

- *nominal volume*
- *50 % of the nominal volume*
- *10 % of the nominal volume*

## In house verification (COVID-19 emergency)

- One point (nominal volume)

## Method: Gravimetric

- A good analytical balance is needed

<sup>1</sup>ISO 8655 - 1 to 7:2002 – *Piston-operated volumetric apparatus*





# Principle of Gravimetric Method

***A COVID-19 EMERGENCY APPROACH***

***“The Mass of water delivered by the pipette is measured”***







# Equipment

## Recommended Balance characteristics \*

<b>Nominal Delivered Volume</b>	<b>Resolution (mg)</b>	<b>Repeatability (mg)</b>	<b>Standard uncertainty (mg)</b>
$1 \mu\text{L} < V < 10 \mu\text{L}$	0,001	0,002	0,002
$10 \mu\text{L} < V \leq 100 \mu\text{L}$	0,01	0,02	0,02
$100 \mu\text{L} < V \leq 10 \text{ mL}$	0,1	0,2	0,2
$10 \text{ mL} < V \leq 200 \text{ mL}$	1	2	2

**\* To be considered necessary**



## Equipment

### **Thermometer 1**

Measurement of the liquid temperature in order to determine its density.

### **Thermometer 2**

Measurement of the air temperature

### **Barometer \***

Measurement of the atmospheric pressure

### **Weighing vessel**

For the weighing of the water mass preferably with a lid to minimize evaporation

\* *If not available, ask for the current local weather report*



# General Techniques

## **Calibration Liquid**

The calibration liquid is normally distilled or de-ionized water

## **Temperature stabilisation**


It is recommended that the piston apparatus and the liquid used in calibration are placed near the balance for 2 hours before calibration

## **Ambient conditions**

It is recommended to control the ambient conditions mainly be as close as possible to 20 °C and have a humidity larger than 50%.



# General Procedure



$$V_{20} = (I_L - I_E) \times (\rho_W - \rho_A)^{-1} \times \left(1 - \frac{\rho_A}{\rho_B}\right) \times [1 - \gamma(t - 20)]$$



## Volume Calculation [1]

$$V_{20} = (I_L - I_E) \times (\rho_W - \rho_A)^{-1} \times \left(1 - \frac{\rho_A}{\rho_B}\right) \times [1 - \gamma(t - 20)]$$

### Where:

- $V_{20}$  = Volume, at reference temperature equal to 20 °C, in ml  
 $I_L$  = Balance reading of vessel with water, in g  
 $I_E$  = Balance reading of empty vessel, in g  
 $\rho_W$  = Density of water at calibration temperature t in g / ml  
 $\rho_A$  = Density of air in g / ml  
 $\rho_B$  = Density of balance weights (8,0 g / ml)  
 $\gamma$  = Coefficient of cubic thermal expansion of the material of which the piston operating apparatus is made, 0,00024 /°C  
 $t$  = Temperature of water used for calibration in °C

[1] ISO 4787: 2010 (E)





## Volume Calculation – Short version

$$V_{20} = (I_L - I_E) \times Z \quad (\text{Eq. 1})$$

$$Z = (\rho_W - \rho_A)^{-1} \times \left( 1 - \frac{\rho_A}{\rho_B} \right) \times [1 - \gamma(t - 20)]$$

To facilitate an easy calculation of the pipette's volume at reference temperature,  $V_{20}$ , using Eq. 1, **the factor Z** has been calculated for polypropylene at common air pressure versus temperature (**Table 1**)



# Volume Calculation – Z factor

Table 1.

**Factor Z for Polypropylene ( $\gamma = 240 \times 10^{-6}$ )**

T [°C]	Water density	Air Density (for relative humidity 50%)											
		Air Pressure [mbar]											
	[kg/m <sup>3</sup> ]	930	940	950	960	970	980	990	1000	1010	1020	1030	1040
15	999,100	1,00309	1,00310	1,00311	1,00312	1,00313	1,00314	1,00315	1,00316	1,00317	1,00318	1,00319	1,00320
15,5	999,023	1,00304	1,00305	1,00306	1,00307	1,00308	1,00310	1,00311	1,00312	1,00313	1,00314	1,00315	1,00316
16	998,943	1,00300	1,00301	1,00302	1,00303	1,00304	1,00305	1,00306	1,00307	1,00309	1,00310	1,00311	1,00312
16,5	998,860	1,00296	1,00297	1,00298	1,00299	1,00300	1,00301	1,00303	1,00304	1,00305	1,00306	1,00307	1,00308
17	998,775	1,00293	1,00294	1,00295	1,00296	1,00297	1,00298	1,00299	1,00300	1,00301	1,00302	1,00303	1,00304
17,5	998,686	1,00289	1,00290	1,00291	1,00292	1,00293	1,00294	1,00296	1,00297	1,00298	1,00299	1,00300	1,00301
18	998,595	1,00286	1,00287	1,00288	1,00289	1,00290	1,00291	1,00292	1,00294	1,00295	1,00296	1,00297	1,00298
18,5	998,502	1,00283	1,00284	1,00285	1,00287	1,00288	1,00289	1,00290	1,00291	1,00292	1,00293	1,00294	1,00295
19	998,405	1,00281	1,00282	1,00283	1,00284	1,00285	1,00286	1,00287	1,00288	1,00289	1,00290	1,00291	1,00292
19,5	998,306	1,00279	1,00280	1,00281	1,00282	1,00283	1,00284	1,00285	1,00286	1,00287	1,00288	1,00289	1,00290
20	998,204	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282	1,00283	1,00284	1,00285	1,00286	1,00287	1,00288
20,5	998,099	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282	1,00283	1,00284	1,00285	1,00286
21	997,992	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282	1,00283	1,00284	1,00285
21,5	997,882	1,00272	1,00273	1,00274	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282	1,00283	1,00284
22	997,770	1,00271	1,00273	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282	1,00283
22,5	997,655	1,00271	1,00272	1,00273	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282
23	997,538	1,00270	1,00271	1,00272	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282
23,5	997,418	1,00270	1,00271	1,00272	1,00273	1,00274	1,00275	1,00276	1,00277	1,00279	1,00280	1,00281	1,00282
24	997,296	1,00270	1,00271	1,00272	1,00273	1,00274	1,00275	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282
24,5	997,171	1,00271	1,00272	1,00273	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00281	1,00282
25	997,044	1,00271	1,00272	1,00273	1,00274	1,00275	1,00276	1,00277	1,00278	1,00279	1,00280	1,00282	1,00283



## In house verification of micropipettes in 9 simple steps \*

1. *Operate the pipette as usually*
2. *Set the pipette volume to be verified*
3. *Measure the water temperature in the water reservoir*
4. *Put the water receiving vessel on the balance and zero the balance*
5. *Aspirate the set volume of water*
6. *Transfer the water to the receiving vessel*
7. *Record balance indication*
8. *Record environmental conditions (Barometric pressure)*
9. ***Use Equation 1 and Table 1 to calculate the volume***



\* *But keep in mind Good Laboratory Practices (see Flyer in last slide)*



## Analysis of results

The obtained error can be compared with the systematic error of the pipette described in ISO 8655-2 or with the previous error described in the calibration certificate.

Calculate the systematic error  $e_S$  of the piston pipette using the formula:

$$e_S = V - \bar{V}$$

Where

$e_S$  is the absolute systematic error, expressed in units of volume;  
 $V_S$  is the selected test volume at the piston pipette under test;  
 $\bar{V}$  is the average result of the volume measurements.



IPO



The European National Metrology Institutes can assist hospital laboratories to identify errors in daily practices, raise awareness and develop best practices and methodologies to use in challenging conditions in order to ensure the reliability of the performed diagnostic tests.

[Watch the video 'Internal verification of piston pipettes'](#)

## METROLOGY FOR MEDICINE & SOCIETY

COVID-19 Crisis & The use of pipettes in diagnostic tests  
THE ROLE OF METROLOGY


### Challenge

The unprecedented conditions public health institutions experience due to COVID-19 pandemic crisis have forced hospital laboratories to perform an enormous number of diagnostic tests daily using pipettes for the handling of the samples.

Correct measurement of sample volume depends on a combination of

- Pipette (Calibration) and
- Operator (Training - GLP)

Incorrect handling by the operator



Operational status of the pipette not verified by calibration

Diagnostic test result may not be reliable

**Tracing the Infection**  
*Probability of false positive vs. false negative status*


Hospital laboratories within COVID-19 crisis are forced to perform diagnostic tests under extraordinary, non-ideal conditions. European Metrology Institutes can assist them to identify errors in daily practice, raise awareness and develop best practices and methodologies to use in challenging conditions in order to ensure the reliability of the performed diagnostic tests.

### How can metrology help?


**Recommendations**

**GOOD LABORATORY PRACTICE (GLP)**


1. Set the volume correctly




2. Attach the tip firmly



3. Pre-wet the tip
4. Work at temperature equilibrium
5. Examine the tip before and after dispensing the sample
6. Use standard mode pipetting
7. Pause consistently after aspiration
8. Pull the pipette straight out



9. Immerse the tip to proper depth



10. Use correct pipette tip
11. Use consistent plunger pressure and speed

**PIPETTE CALIBRATION**

1. Frequent calibration
2. Accuracy verification
3. Compliance with specifications
4. Be aware of wear due to excessive use



NATIONAL QUALITY INFRASTRUCTURE SYSTEM  
HELLENIC INSTITUTE OF METROLOGY



Contact  
Ela Bektas (IPO, Portugal) eabektas@ipq.pt  
Zoe Melarobou (EIM, Greece) zoe@eim.gr





## Copyright information

*The copyright for the this material may be held by the authors .The material is provided for reading and for guidance in your work.*

*You may use this material for internal training purposes, as long as due reference is made to the authors and their institutional affiliation.*

*By using this document, you confirm that you have read and agreed with the copyright information.*

*If your organisation would like to issue its own version or translate the material, please contact the copyright holder.*



Thank you!

For more information please  
contact us by e-mail:

[ebatista@ipq.pt](mailto:ebatista@ipq.pt)

[zoe@eim.gr](mailto:zoe@eim.gr)